

# FcRn-Specific antibody



Catalog Number: 110606

## Product name

Polyclonal

FcRn-Specific antibody

## Ig Type

## Specificity

Rabbit IgG

Human; other species not tested.

## Applications

## Antibody description

ELISA, WB

FcRn-Specific Rabbit Polyclonal antibody. Positive WB detected in human placenta tissue, human liver tissue. Observed molecular weight by Western-blot: 35 kDa

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_004107).  
Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



human placenta tissue were subjected to SDS PAGE followed by western blot with Catalog No:110606(FCGRT-Specific antibody) at dilution of 1:300