

FER antibody

Catalog Number: 110619

Product name

FER antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FER Rabbit Polyclonal antibody. Positive FC detected in Jurkat cells. Positive IHC detected in human colon cancer tissue, human colon tissue, human kidney tissue. Positive IF detected in MCF-7 cells. Positive WB detected in NIH/3T3 cells, A431 cells, HeLa cells, MCF-7 cells. Observed molecular weight by Western-blot: 95 kDa

Preparation

This antibody was obtained by immunization of FER recombinant protein (Accession Number: BC017060). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, FC, IF

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IHC: 1:20-1:200

IF: 1:20-1:200

Validations




NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110619(FER Antibody) at dilution of 1:1000




Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:110619(FER Antibody) at dilution of 1:50 (under 10x lens)

FER antibody


Catalog Number: 110619



Immunohistochemistry of paraffin-embedded human colon cancer tissue slide using Catalog No:110619(FER Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of MCF-7 cells using Catalog No:110619(FER Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



1×10^6 Jurkat cells were stained with 0.5ug FER antibody (Catalog No:110619, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.