

FEV antibody

Catalog Number: 110625

Product name

FEV antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FEV Rabbit Polyclonal antibody. Positive WB detected in fetal human brain tissue, mouse brain tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of FEV recombinant protein (Accession Number: NM_017521). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



fetal human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110625(FEV Antibody) at dilution of 1:1000



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110625(FEV Antibody) at dilution of 1:1000