

FFAR2 antibody

Catalog Number: 110626

Product name

FFAR2 antibody

Specificity

Human; other species not tested.

Antibody description

FFAR2 Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue, human colon tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005306).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:110626(FFAR2 Antibody) at dilution of 1:100 (under 10x lens).



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:110626(FFAR2 Antibody) at dilution of 1:100 (under 40x lens).