

FGF18 antibody



Catalog Number: 110634

Product name

FGF18 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FGF18 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, mouse heart tissue. Positive IF detected in HepG2 cells.

Preparation

This antibody was obtained by immunization of FGF18 recombinant protein (Accession Number: NM_033649). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110634(FGF18 antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:110634(FGF18 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)