

# FGF5 antibody



Catalog Number: 110635

## Product name

FGF5 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

## Specificity

Human, Mouse; other species not tested.

IHC: 1:50-1:500

IF: 1:20-1:200

## Antibody description

FGF5 Rabbit Polyclonal antibody. Positive IP detected in A375 cells. Positive WB detected in A375 cells, NIH/3T3 cells, SH-SY5Y cells. Positive IHC detected in human stomach tissue, human placenta tissue. Positive IF detected in NIH/3T3 cells. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of FGF5 recombinant protein (Accession Number: NM\_004464). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IF, IP, IHC


## Dilutions

Recommended Dilution:

## Validations




A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:110635(FGF5 Antibody) at dilution of 1:600




Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:110635(FGF5 Antibody) at dilution of 1:200 (under 10x lens).

# FGF5 antibody


Catalog Number: 110635



Immunohistochemistry of paraffin-embedded human stomach tissue slide using Catalog No:110635(FGF5 Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of NIH/3T3 cells using Catalog No:110635(FGF5 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



IP Result of anti-FGF5 (IP:Catalog No:110635, 4ug; Detection:Catalog No:110635 1:500) with A375 cells lysate 3200ug.