

FGFBP2 antibody



Catalog Number: 110638

Product name

FGFBP2 antibody

Specificity

Human; other species not tested.

Antibody description

FGFBP2 Rabbit Polyclonal antibody. Positive IF detected in BxPC-3 cells. Positive WB detected in SH-SY5Y cells, BxPC-3 cells. Observed molecular weight by Western-blot: 37kd

Preparation

This antibody was obtained by immunization of FGFBP2 recombinant protein (Accession Number: NM_031950). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IF: 1:10-1:100

Validations



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with Catalog No:110638(FGFBP2 antibody) at dilution of 1:400



Immunofluorescent analysis of BxPC-3 cells using Catalog No:110638(FGFBP2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)