

# FGFR1OP antibody



Catalog Number: 110640

## Product name

FGFR1OP antibody

## Immunogen

[Human FGFR1OP Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FGFR1OP Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue. Positive IHC detected in human lymphoma tissue. Observed molecular weight by Western-blot: 50kd

## Preparation

This antibody was obtained by immunization of FGFR1OP recombinant protein (Accession Number: NM\_194429). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110640(FGFR1OP antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:110640(FGFR1OP antibody) at dilution of 1:50 (under 10x lens)