

FGFR2 antibody

Catalog Number: 110642

Product name

FGFR2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FGFR2 Rabbit Polyclonal antibody. Positive WB detected in A549 cells, A431 cells, HeLa cells, HepG2 cells, SGC-7901 cells. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 120 kDa, 85-92 kDa

Preparation

This antibody was obtained by immunization of FGFR2 recombinant protein (Accession Number: NM_001144918). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:110642(FGFR2 antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells using Catalog No:110642(FGFR2 Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)