

FGL2 antibody

Catalog Number: 110645

Product name

FGL2 antibody

IHC: 1:20-1:200

IF: 1:10-1:100


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

FGL2 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive IF detected in HepG2 cells. Positive WB detected in human ileum tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 70 kDa, 50-60 kDa




Preparation

This antibody was obtained by immunization of FGL2 recombinant protein (Accession Number: BC017813). Purification method: Antigen affinity purified.

human ileum tissue were subjected to SDS PAGE followed by western blot with Catalog No:110645(FGL2 antibody) at dilution of 1:400

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110645(FGL2 antibody) at dilution of 1:50 (under 10x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IF



Dilutions

Recommended Dilution:


WB: 1:500-1:5000

Immunohistochemical of paraffin-embedded

FGL2 antibody

Catalog Number: 110645

human colon cancer using Catalog No:110645(FGL2 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using FGL2 antibody Catalog No:110645 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red). Blue pseudocolor = DAPI (fluorescent DNA dye).