

Fhl1 antibody



Catalog Number: 110648

Product name

Fhl1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Fhl1 Rabbit Polyclonal antibody. Positive WB detected in human heart tissue. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_010211). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110648(Fhl1 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:110648(Fhl1 antibody) at dilution of 1:50 (under 10x lens)