

# FHL1 antibody

Catalog Number: 110649

## Product name

FHL1 antibody

## Immunogen

[Human FHL1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Zebrafish; other species not tested.

## Antibody description

FHL1 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, mouse colon tissue, mouse skeletal muscle tissue. Positive IP detected in mouse skeletal muscle tissue. Positive IHC detected in human ovary tumor tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of FHL1 recombinant protein (Accession Number: NM\_001159700). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IP

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000


IP: 1:200-1:2000

IHC: 1:20-1:200

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110649(FHL1 antibody) at dilution of 1:1000




Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:110649(FHL1 antibody) at dilution of 1:100 (under 10x lens)

# FHL1 antibody

Catalog Number: 110649

---



IP Result of anti-FHL1 (IP:Catalog No:110649, 3ug; Detection:Catalog No:110649 1:700) with mouse skeletal muscle tissue lysate 3600ug.