

FHL5 antibody



Catalog Number: 110652

Product name

FHL5 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FHL5 Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, human heart tissue. Observed molecular weight by Western-blot: 33 kDa

Preparation

This antibody was obtained by immunization of FHL5 recombinant protein (Accession Number: XM_047419567). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110652(FHL5 Antibody) at dilution of 1:2000