

FKBP10 antibody



Catalog Number: 110669

Product name

FKBP10 antibody

Recommended Dilution:

WB: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IP: 1:200-1:2000

IHC: 1:20-1:200

Antibody description

FKBP10 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive IF detected in HEK-293 cells. Positive IP detected in HEK-293 cells. Positive WB detected in HEK-293 cells, fetal human brain tissue, human brain tissue, human kidney tissue, human placenta tissue, NIH/3T3 cells. Observed molecular weight by Western-blot: 64 kDa

IF: 1:50-1:500

Preparation

This antibody was obtained by immunization of FKBP10 recombinant protein (Accession Number: NM_021939). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications


ELISA, IHC, WB, IP, IF

Dilutions

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110669(FKBP10 Antibody) at dilution of 1:600




Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:110669(FKBP10 antibody) at dilution of 1:50 (under 10x lens)


FKBP10 antibody



Catalog Number: 110669



Anti-Rabbit IgG(H+L)



Immunofluorescent analysis of (10% Formaldehyde) fixed HEK-293 cells using Catalog No:110669(FKBP10 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat

IP Result of anti-FKBP10 (IP:Catalog No:110669, 4ug; Detection:Catalog No:110669 1:600) with HEK-293 cells lysate 2800ug.