

# FKBP14 antibody



Catalog Number: 110670

---

## Product name

FKBP14 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FKBP14 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, human brain tissue. Observed molecular weight by Western-blot: 28 kDa

## Preparation

This antibody was obtained by immunization of FKBP14 recombinant protein (Accession Number: NM\_017946). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110670(FKBP14 antibody) at dilution of 1:500