

# FKBP3 antibody



Catalog Number: 110674

---

## Product name

FKBP3 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

FKBP3 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, A431 cells, human brain tissue, PC-3 cells. Observed molecular weight by Western-blot: 25kd

## Preparation

This antibody was obtained by immunization of FKBP3 recombinant protein (Accession Number: NM\_002013). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110674(FKBP3 antibody) at dilution of 1:500