

FLNB antibody



Catalog Number: 110690

Product name

FLNB antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FLNB Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HeLa cells, NIH/3T3 cells. Observed molecular weight by Western-blot: 250-278kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001457).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol
pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110690(FLNB antibody) at dilution of 1:300



Immunofluorescent analysis of HepG2 cells, using FLNB antibody Catalog No:110690 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).