

# FOXB1 antibody



Catalog Number: 110691

## Product name

FOXB1 antibody

## Immunogen

[Human FOXB1 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FOXB1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, mouse liver tissue. Observed molecular weight by Western-blot: 35 kDa

## Preparation

This antibody was obtained by immunization of FOXB1 recombinant protein (Accession Number: NM\_012182). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110691(FOXB1



Antibody) at dilution of 1:300

mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110691(FOXB1 Antibody) at dilution of 1:300