

FOXD4L6 antibody



Catalog Number: 110695

Product name

FOXD4L6 antibody

Specificity

Human; other species not tested.

Antibody description

FOXD4L6 Rabbit Polyclonal antibody. Positive WB detected in A431 cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 46 kDa

Preparation

This antibody was obtained by immunization of FOXD4L6 recombinant protein (Accession Number: NM_001085476). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:110695(FOXD4L6 antibody) at dilution of 1:1000