

FOXG1 antibody



Catalog Number: 110697

Product name

FOXG1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FOXG1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, mouse brain tissue, SH-SY5Y cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of FOXG1 recombinant protein (Accession Number: BC035020). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110697(FOXG1 antibody) at dilution of 1:500