

FOXJ3 antibody



Catalog Number: 110699

Product name

FOXJ3 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

FOXJ3 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, HepG2 cells. Observed molecular weight by Western-blot: 66-69kd

Preparation

This antibody was obtained by immunization of FOXJ3 recombinant protein (Accession Number: XM_047449478). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:110699(FOXJ3 antibody) at dilution of 1:500