

# FOXL2-Specific antibody



Catalog Number: 110701

## Product name

FOXL2-Specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FOXL2-Specific Rabbit Polyclonal antibody.  
Positive WB detected in mouse ovary tissue.  
Observed molecular weight by Western-blot: 49 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_023067).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:110701(FOXL2-Specific Antibody) at dilution of 1:300