

# FOXM1 antibody

Catalog Number: 110702

## Product name

FOXM1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FOXM1 Rabbit Polyclonal antibody. Positive WB detected in mouse thymus tissue, human testis tissue, L02 cells, mouse colon tissue, mouse spleen tissue, mouse testis tissue. Positive IP detected in L02 cells. Observed molecular weight by Western-blot: 110kd

## Preparation

This antibody was obtained by immunization of FOXM1 recombinant protein (Accession Number: NM\_202003). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



mouse thymus tissue were subjected to SDS PAGE followed by western blot with Catalog No:110702(FOX M1 antibody) at dilution of 1:1000



IP Result of anti-FOX M1 (IP:Catalog No:110702, 5ug; Detection:Catalog No:110702 1:700) with L02 cells lysate 1560ug.