

FOXO4 antibody



Catalog Number: 110705

Product name

FOXO4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FOXO4 Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in human heart tissue, mouse heart tissue. Observed molecular weight by Western-blot: 45-48kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005938).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:10-1:100

Validations



human heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110705(FOXO4 antibody) at dilution of 1:800



Immunofluorescent analysis of HepG2 cells, using FOXO4 antibody Catalog No:110705 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).