

FOXO6 antibody

Catalog Number: 110707

Product name

FOXO6 antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FOXO6 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Positive IHC detected in human gliomas tissue, human brain tissue. Observed molecular weight by Western-blot: 70 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: AI361654). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110707(FOXO6 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110707(FOXO6 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:110707(FOXO6

FOXO6 antibody

Catalog Number: 110707

antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-FOXO6 (IP:Catalog No:110707, 3ug; Detection:Catalog No:110707 1:300) with mouse brain tissue lysate 4000ug.

