

# FOXP2 antibody

Catalog Number: 110709

## Product name

FOXP2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FOXP2 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse brain tissue, mouse spleen tissue. Positive IP detected in mouse spleen tissue. Observed molecular weight by Western-blot: 70-80 kDa

## Preparation

This antibody was obtained by immunization of FOXP2 recombinant protein (Accession Number: BC018016). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110709(FOXP2 antibody) at dilution of 1:300



IP Result of anti-FOXP2 (IP:Catalog No:110709, 4ug; Detection:Catalog No:110709 1:500) with mouse spleen tissue lysate 4000ug.