

FOXP3 antibody



Catalog Number: 110710

Product name

FOXP3 antibody

Immunogen

[Human FOXP3 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

FOXP3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A375 cells, HeLa cells, RAW264.7 cells. Observed molecular weight by Western-blot: 47-50 kDa

Preparation

This antibody was obtained by immunization of FOXP3 recombinant protein (Accession Number: NM_014009). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:200-1:2000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110710(FOXP3 antibody) at dilution of 1:1000