

FOXQ1 antibody



Catalog Number: 110713

Product name

FOXQ1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FOXQ1 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, rat kidney tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of FOXQ1 recombinant protein (Accession Number: NM_033260). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110713(FOXQ1 antibody) at dilution of 1:800