

Frizzled 10 antibody



Catalog Number: 110721

Product name

Frizzled 10 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Frizzled 10 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HeLa cells. Positive IHC detected in human kidney tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of Frizzled 10 recombinant protein (Accession Number: NM_007197). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110721(FZD10 antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human kidney using Catalog No:110721(FZD10 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:110721(FZD10 antibody) at dilution of 1:100 (under 40x lens)