

Flotillin 2 antibody



Catalog Number: 110723

Product name

Flotillin 2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

Flotillin 2 Rabbit Polyclonal antibody. Positive WB detected in A375 cells, rat lung tissue. Observed molecular weight by Western-blot: 42 kDa

Preparation

This antibody was obtained by immunization of Flotillin 2 recombinant protein (Accession Number: BC017292). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



A375 cells were subjected to SDS PAGE followed by western blot with Catalog No:110723(FLOT2 antibody) at dilution of 1:300