

FLT3 antibody

Catalog Number: 110724

Product name

FLT3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FLT3 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 130-140 kDa

Preparation

This antibody was obtained by immunization of FLT3 recombinant protein (Accession Number: BC036028). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:110724(FLT3 antibody) at dilution of 1:500



IP Result of anti-FLT3 (IP:Catalog No:110724, 4ug; Detection:Catalog No:110724 1:500) with mouse brain tissue lysate 2640ug.