

FMO5 antibody



Catalog Number: 110730

Product name

FMO5 antibody

Immunogen

[Human FMO5 Recombinant protein \(His tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FMO5 Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, mouse lung tissue. Positive IP detected in mouse lung tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of FMO5 recombinant protein (Accession Number: NM_001144830). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110730(FMO5 antibody) at dilution of 1:1200



IP Result of anti-FMO5 (IP:Catalog No:110730, 4ug; Detection:Catalog No:110730 1:500) with mouse lung tissue lysate 4000ug.