

# FMR1NB antibody



Catalog Number: 110733

## Product name

Polyclonal

FMR1NB antibody

## Ig Type

## Specificity

Rabbit IgG

human; other species not tested.

## Applications

## Antibody description

ELISA, WB

FMR1NB Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, A549 cells, SKOV-3 cells. Observed molecular weight by Western-blot: 66kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of FMR1NB recombinant protein (Accession Number: NM\_152578). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110733(FMR1NB antibody) at dilution of 1:1000