

FNDC3A antibody



Catalog Number: 110736

Product name

FNDC3A antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FNDC3A Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, mouse placenta tissue. Observed molecular weight by Western-blot: 120-130 kDa

Preparation

This antibody was obtained by immunization of FNDC3A recombinant protein (Accession Number: XM_017020440). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110736(FNDC3A Antibody) at dilution of 1:1000