

FNDC3B antibody

Catalog Number: 110737

Product name

FNDC3B antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FNDC3B Rabbit Polyclonal antibody. Positive IP detected in HepG2 cells. Positive WB detected in Jurkat cells, HEK-293 cells, HepG2 cells, HL-60 cells. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 150 kDa, 70 kDa

Preparation

This antibody was obtained by immunization of FNDC3B recombinant protein (Accession Number: XM_024453716). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IP, IHC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:110737(FNDC3B antibody) at dilution of 1:2000



IP Result of anti-FNDC3B (IP:Catalog No:110737, 4ug; Detection:Catalog No:110737 1:1000) with HepG2 cells lysate 2400ug.




Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog

FNDC3B antibody

Catalog Number: 110737

No:110737(FNDC3B Antibody) at dilution of 1:200
(under 10x lens).



Immunohistochemistry of paraffin-embedded
human liver tissue slide using Catalog
No:110737(FNDC3B Antibody) at dilution of 1:200
(under 40x lens).