

FNDC5 antibody

Catalog Number: 110738

Product name

FNDC5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

FNDC5 Rabbit Polyclonal antibody. Positive IHC detected in human skeletal muscle tissue. Positive IP detected in mouse liver tissue. Positive WB detected in mouse placenta tissue, mouse brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 25kDa - 30 kDa

Preparation

This antibody was obtained by immunization of FNDC5 recombinant protein (Accession Number: NM_001171941). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

IP: 1:200-1:1000

IHC: 1:20-1:200

Validations



mouse liver, mouse skeletal muscle and rat liver were subjected to SDS PAGE followed by western blot with Catalog No:110738(FNDC5 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours



Immunohistochemistry of paraffin-embedded human skeletal muscle slide using Catalog No:110738(FNDC5 Antibody) at dilution of 1:50