

# FXYD2 antibody



Catalog Number: 110754

## Product name

FXYD2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FXYD2 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in HEK-293 cells, human kidney tissue, human liver tissue, human skeletal muscle tissue, mouse skeletal muscle tissue. Observed molecular weight by Western-blot: 7 kDa

## Preparation

This antibody was obtained by immunization of FXYD2 recombinant protein (Accession Number: NM\_001680). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110754(FXYD2 antibody) at dilution of 1:1000



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:110754(FXYD2 antibody) at dilution of 1:100 (under 10x lens)