

# FZR1 antibody



Catalog Number: 110761

## Product name

FZR1 antibody

WB: 1:200-1:2000

IP: 1:200-1:2000

## Specificity

Human, Mouse, Rat; other species not tested.

IF: 1:50-1:500

## Antibody description

FZR1 Rabbit Polyclonal antibody. Positive IP detected in mouse heart tissue. Positive WB detected in NIH/3T3 cells, C6 cells, HeLa cells, HepG2 cells, Jurkat cells, mouse heart tissue, mouse liver tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of FZR1 recombinant protein (Accession Number: NM\_016263). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB, IP, IF


## Dilutions

Recommended Dilution:


## Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:110761(FZR1 Antibody) at dilution of 1:600



IP Result of anti-FZR1 (IP:Catalog No:110761, 4ug; Detection:Catalog No:110761 1:500) with mouse heart tissue lysate 3200ug.



# FZR1 antibody



Catalog Number: 110761

---

Immunofluorescent analysis of (10% Formaldehyde) fixed HeLa cells using Catalog No:110761(FZR1 Antibody) at dilution of 1:50 and

Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)