

G0S2 antibody

Catalog Number: 110762

Product name

G0S2 antibody

Immunogen

[Human G0S2 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

G0S2 Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells. Positive IF detected in mouse (embryonic day 18) skin tissue.

Preparation

This antibody was obtained by immunization of G0S2 recombinant protein (Accession Number: NM_015714). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:25 - 1:100

Validations



IF result of G0S2 antibody (Catalog No:110762,1:100, Green) with skin samples (mouse) at embryonic day 18 by Dr. Colin Jahoda.



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110762(G0S2 Antibody) at dilution of 1:1000