

# G6PC3 antibody



Catalog Number: 110765

---

## Product name

G6PC3 antibody

## Specificity

Human; other species not tested.

## Antibody description

G6PC3 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue.

## Preparation

This antibody was obtained by immunization of G6PC3 recombinant protein (Accession Number: NM\_138387). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:110765(G6PC3 Antibody) at dilution of 1:200 (under 40x lens).