

G6PC3 antibody



Catalog Number: 110765

Product name

G6PC3 antibody

Specificity

Human; other species not tested.

Antibody description

G6PC3 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue.

Preparation

This antibody was obtained by immunization of G6PC3 recombinant protein (Accession Number: NM_138387). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human liver tissue slide using Catalog No:110765(G6PC3 Antibody) at dilution of 1:200 (under 40x lens).