

GAB2 antibody

Catalog Number: 110768

Product name

GAB2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GAB2 Rabbit Polyclonal antibody. Positive IHC detected in human breast cancer tissue, human tonsillitis tissue. Positive WB detected in MCF-7 cells, mouse spleen tissue. Observed molecular weight by Western-blot: 70-74 kDa

Preparation

This antibody was obtained by immunization of GAB2 recombinant protein (Accession Number: NM_080491). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



MCF-7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110768(GAB2 Antibody) at dilution of 1:600



Immunohistochemistry of paraffin-embedded human breast cancer tissue slide using Catalog No:110768(GAB2 Antibody) at dilution of 1:50 (under 40x lens)