

# Frizzled 8 antibody

Catalog Number: 110774

## Product name

Frizzled 8 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

Frizzled 8 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse embryo tissue, mouse heart tissue, mouse lung tissue. Observed molecular weight by Western-blot: 70-73 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_031866). Purification method: Antigen affinity purified.

## Formulation

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions


Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110774(FZD8 Antibody) at dilution of 1:600



WB result of Catalog No:110774 (FZD8 Antibody).