

# FRK antibody

Catalog Number: 110776

## Product name

FRK antibody

Recommended Dilution:

WB: 1:500-1:5000

## Specificity

Human, Mouse, Rat; other species not tested.


IHC: 1:20-1:200

IF: 1:10-1:100

## Antibody description

FRK Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, K-562 cells, MCF7 cells, mouse kidney tissue. Positive IHC detected in human breast cancer tissue, human bladder tissue, human kidney tissue. Positive IF detected in HepG2 cells. Positive FC detected in HepG2 cells. Observed molecular weight by Western-blot: 54-57 kDa

## Validations



## Preparation

This antibody was obtained by immunization of FRK recombinant protein (Accession Number: XM\_005266880). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110776(FRK antibody) at dilution of 1:800

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal


## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, IF, FC


## Dilutions




Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110776(FRK antibody) at dilution of 1:50 (under 10x lens)

# FRK antibody


Catalog Number: 110776



Immunohistochemical of paraffin-embedded human breast cancer using Catalog No:110776(FRK antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using FRK antibody Catalog No:110776 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).



1X10<sup>6</sup> HepG2 cells were stained with 0.2ug FRK antibody (Catalog No:110776, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.