

FUBP3 antibody



Catalog Number: 110789

Product name

FUBP3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FUBP3 Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in HEK-293 cells.

Preparation

This antibody was obtained by immunization of FUBP3 recombinant protein (Accession Number: BC007874). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human prostate cancer using Catalog No:110789(FUBP3 antibody) at dilution of 1:100



(under 25x lens)

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110789(FUBP3 Antibody) at dilution of 1:600