

FUK antibody



Catalog Number: 110792

Product name

FUK antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

FUK Rabbit Polyclonal antibody. Positive IHC detected in mouse kidney tissue. Positive WB detected in mouse kidney tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 107kd

Preparation

This antibody was obtained by immunization of FUK recombinant protein (Accession Number: XM_011522929). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110792(FUK antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded mouse kidney tissue slide using Catalog No:110792(FUK Antibody) at dilution of 1:50 (under 40x lens)