

FURIN antibody

Catalog Number: 110794

Product name

FURIN antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


Specificity

Human, Mouse, Rat; other species not tested.

Validations

Antibody description

FURIN Rabbit Polyclonal antibody. Positive IHC detected in human placenta tissue. Positive WB detected in HeLa cells, COLO 320 cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 90 kDa




Preparation

This antibody was obtained by immunization of FURIN recombinant protein (Accession Number: NM_001382619). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110794(FURIN Antibody) at dilution of 1:600

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

IP Result of anti-FURIN (IP:Catalog No:110794, 4ug; Detection:Catalog No:110794 1:500) with HeLa cells lysate 1400ug.

Ig Type

Rabbit IgG


Applications

ELISA, WB, IHC, IP

Dilutions

Recommended Dilution:

WB: 1:200-1:2000




Immunohistochemistry of paraffin-embedded human placenta tissue slide using Catalog

FURIN antibody

Catalog Number: 110794

No:110794(FURIN Antibody) at dilution of 1:200
(under 10x lens).



Immunohistochemistry of paraffin-embedded
human placenta tissue slide using Catalog
No:110794(FURIN Antibody) at dilution of 1:200
(under 40x lens).