

# FUSIP1 antibody



Catalog Number: 110796

## Product name

Polyclonal

FUSIP1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

FUSIP1 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 31 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of FUSIP1 recombinant protein (Accession Number: NM\_006625). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110796(FUSIP1 Antibody) at dilution of 1:1000