

FUT4 antibody

Catalog Number: 110799

Product name

FUT4 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human; other species not tested.

Antibody description

FUT4 Rabbit Polyclonal antibody. Positive IHC detected in human liver cancer tissue, human gliomas tissue, human skin cancer tissue. Positive WB detected in HeLa cells, HepG2 cells, HL-60 cells, Jurkat cells. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 95-110 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002033). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, IP, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:2000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:110799(FUT4 antibody) at dilution of 1:500



IP Result of anti-FUT4 (IP:Catalog No:110799, 4ug; Detection:Catalog No:110799 1:500) with HeLa cells lysate 2480ug.




Immunohistochemistry of paraffin-embedded human liver cancer tissue slide using Catalog

FUT4 antibody

Catalog Number: 110799

No:110799(FUT4 Antibody) at dilution of 1:100
(under 10x lens)



Immunohistochemistry of paraffin-embedded
human liver cancer tissue slide using Catalog
No:110799(FUT4 Antibody) at dilution of 1:100
(under 40x lens)