

# GABRB1 antibody



Catalog Number: 110810

## Product name

GABRB1 antibody

## Immunogen

[Human GABRB1 Recombinant protein \(GST tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

GABRB1 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Positive IP detected in mouse testis tissue. Observed molecular weight by Western-blot: 50-54 kDa

## Preparation

This antibody was obtained by immunization of GABRB1 recombinant protein (Accession Number: NM\_000812). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:110810(GABRB1 antibody) at dilution of 1:800



IP Result of anti-GABRB1 (IP:Catalog No:110810, 4ug; Detection:Catalog No:110810 1:500) with mouse testis tissue lysate 4000ug.